


Silimed Smooth Round Silicone Gel Implants



Silimed Smooth Round Silicone Gel Implants


The Silimed® Smooth Round Silicone Gel Implants are comprised of a silicone-elastomer shell, smooth-surfaced with high strength silicone gel. Provided in two profiles: HIGH and MODERATE. Supplied sterile.

HIGH PROFILE



REF	VOLUME cc	DIMENSIONS cm	
		A	B
10521-095	95	8.2	2.5
10521-115	115	8.8	2.6
10521-135	135	9.2	2.7
10521-155	155	9.7	2.9
10521-175	175	9.9	3.0
10521-195	195	10.1	3.1
10521-215	215	10.2	3.4
10521-235	235	10.7	3.4
10521-255	255	11.0	3.5
10521-275	275	11.1	3.7
10521-295	295	11.3	3.8
10521-315	315	11.4	4.0
10521-335	335	11.8	4.1
10521-355	355	11.9	4.2
10521-385	385	12.3	4.3
10521-415	415	12.7	4.4
10521-455	455	12.8	4.5
10521-495	495	13.3	4.7
10521-545	545	13.9	4.8
10521-605	605	14.3	4.8
10521-655	655	15.1	4.8
10521-695	695	15.4	4.9
10521-745	745	14.8	6.1
10521-795	795	15.1	6.2
10521-845	845	15.3	6.3
10521-895	895	15.8	6.4
10521-945	945	16.1	6.5
10521-995	995	16.4	6.6

MODERATE PROFILE



REF	VOLUME cc	DIMENSIONS cm	
		A	B
10512-080	80	8.1	2.1
10512-100	100	8.6	2.3
10512-120	120	9.2	2.4
10512-140	140	9.5	2.6
10512-160	160	10.0	2.7
10512-180	180	10.5	2.8
10512-200	200	10.8	2.9
10512-220	220	11.3	2.9
10512-240	240	11.6	3.0
10512-260	260	11.9	3.1
10512-280	280	12.2	3.2
10512-300	300	12.5	3.3
10512-330	330	12.7	3.4
10512-360	360	13.0	3.6
10512-390	390	13.2	3.8
10512-420	420	13.6	3.9
10512-450	450	14.0	4.0
10512-500	500	14.3	4.1
10512-550	550	14.8	4.3
10512-600	600	15.4	4.4
10512-650	650	15.8	4.5
10512-700	700	16.1	4.7
10512-750	750	16.0	5.0
10512-800	800	16.1	5.4
10512-850	850	16.2	5.7
10512-900	900	16.4	6.0
10512-950	950	16.6	6.1
10512-1000	1000	17.0	6.2



Caution: Investigational Device limited by Federal (USA) law to investigational use.


Silimed Textured Round Silicone Gel Implants



Silimed Textured Round Silicone Gel Implants


The Silimed® Textured Round Silicone Gel Implants are comprised of a high strength silicone gel-filled implant whose envelope presents a textured surface with open pores. Provided in two profiles: HIGH and LOW. Supplied sterile, with a polyethylene sleeve for insertion.

HIGH PROFILE



REF	VOLUME cc	DIMENSIONS cm	
		A	B
20621-95	95	7.7	2.9
20621-115	115	8.1	3.2
20621-135	135	8.7	3.3
20621-155	155	8.8	3.6
20621-175	175	9.3	3.7
20621-195	195	9.7	3.9
20621-215	215	10.0	3.9
20621-235	235	10.5	3.9
20621-255	255	11.0	3.9
20621-285	285	11.4	4.0
20621-305	305	11.7	4.0
20621-325	325	11.8	4.4
20621-355	355	12.1	4.4
20621-385	385	12.5	4.5
20621-435	435	13.2	4.6
20621-485	485	13.6	4.8
20621-525	525	13.8	5.0
20621-575	575	14.1	5.2
20621-625	625	14.6	5.3
20621-695	695	15.1	5.5
20621-755	755	15.6	5.9
20621-815	815	16.0	6.0
20621-880	880	16.4	6.1
20621-935	935	16.8	6.3
20621-995	995	17.1	6.4

LOW PROFILE



REF	VOLUME cc	DIMENSIONS cm	
		A	B
20610-40	40	7.0	1.8
20610-60	60	7.3	2.1
20610-80	80	8.0	2.2
20610-100	100	8.6	2.3
20610-120	120	9.1	2.4
20610-140	140	10.1	2.4
20610-160	160	10.5	2.5
20610-190	190	11.4	2.6
20610-220	220	12.1	2.7
20610-250	250	12.3	2.8
20610-280	280	13.0	2.9
20610-310	310	13.4	3.0
20610-340	340	13.7	3.1
20610-370	370	14.5	3.1
20610-400	400	14.8	3.2
20610-450	450	15.4	3.2
20610-500	500	15.9	3.3
20610-550	550	16.5	3.5
20610-600	600	17.0	3.7
20610-650	650	17.5	3.8
20610-700	700	17.9	3.8
20610-750	750	18.3	4.0
20610-800	800	18.7	4.1
20610-850	850	19.1	4.2
20610-900	900	19.5	4.3



Caution: Investigational Device limited by Federal (USA) law to investigational use.


Silimed Textured Shaped Silicone Gel Implants - Nuance® & Enhance®



Silimed Textured Shaped Silicone Gel Implants - Nuance® & Enhance®


The Silimed® Nuance® & Enhance® Textured Shaped Silicone Gel Implants are high strength silicone gel implants whose envelopes present a textured surface with open pores, and feature an orientation mark for positioning. The Nuance® Inferior Pole Implant is provided in three profiles: HIGH, MODERATE, and LOW. The Enhance® is a Superior Pole Implant. Supplied sterile, with a polyethylene sleeve for insertion.

NUANCE® HIGH PROFILE




REF	VOLUME cc	DIMENSIONS cm		
		A	B	C
20646-180	180	9.8	8.3	4.3
20646-240	240	10.5	8.9	4.7
20646-300	300	11.4	10.1	5.0
20646-370	370	12.5	10.8	5.3
20646-480	480	13.5	11.5	5.7
20646-550	550	14.6	12.4	5.8

NUANCE® LOW PROFILE




REF	VOLUME cc	DIMENSIONS cm		
		A	B	C
20644-170	170	11.3	9.8	2.8
20644-220	220	12.3	10.6	3.1
20644-270	270	13.3	11.1	3.3
20644-320	320	14.3	12.1	3.5
20644-410	410	15.4	13.1	3.7
20644-500	500	16.3	13.8	4.1

NUANCE® MODERATE PROFILE



REF	VOLUME cc	DIMENSIONS cm		
		A	B	C
20645-120	120	8.9	8.0	3.4
20645-150	150	9.4	8.5	3.6
20645-180	180	10.0	8.9	3.9
20645-210	210	11.0	9.6	4.0
20645-250	250	12.0	10.3	4.2
20645-290	290	12.3	10.7	4.5
20645-330	330	13.0	11.2	4.5
20645-370	370	13.4	11.4	4.7
20645-400	400	13.7	11.8	4.9
20645-450	450	14.1	12.2	5.1
20645-500	500	14.7	12.8	5.3
20645-550	550	15.2	13.1	5.4
20645-600	600	16.0	13.8	5.5
20645-650	650	16.6	14.3	5.9
20645-700	700	16.9	14.5	5.9
20645-750	750	17.3	15.0	6.0
20645-800	800	17.6	15.2	6.2
20645-850	850	17.8	15.5	6.3
20645-900	900	18.3	16.0	6.3
20645-950	950	18.5	16.2	6.4
20645-1000	1000	18.5	16.0	6.4

ENHANCE®



REF	VOLUME cc	DIMENSIONS cm		
		A	B	C
20676-115	115	9.0	8.0	3.2
20676-135	135	9.5	8.5	3.3
20676-155	155	9.9	9.0	3.5
20676-175	175	10.3	9.3	3.6
20676-195	195	10.7	9.8	3.7
20676-225	225	11.1	10.0	4.0
20676-245	245	11.5	10.4	4.1
20676-265	265	11.9	10.9	4.1
20676-285	285	12.3	11.3	4.1
20676-305	305	12.7	11.6	4.1
20676-325	325	13.0	11.9	4.1
20676-365	365	13.4	12.4	4.3
20676-390	390	14.1	13.0	4.3
20676-430	430	14.5	13.4	4.3
20676-480	480	15.2	14.1	4.4
20676-530	530	15.5	14.5	4.5
20676-570	570	16.0	14.8	4.6
20676-640	640	16.2	15.2	4.9
20676-690	690	17.2	16.1	5.0
20676-740	740	17.6	16.4	5.2
20676-790	790	17.9	16.7	5.4
20676-840	840	18.3	17.1	5.5
20676-890	890	18.5	17.4	5.6
20676-940	940	19.0	17.7	5.7
20676-990	990	19.1	18.0	5.8



Caution: Investigational Device limited by Federal (USA) law to investigational use.



PO Box 1490 | Santa Barbara, CA 93116-1490 | Tel: 805.562.3500 | Toll Free: 888.708.0808
Fax: 805.562.8401 | www.sientra.com

MDOC-DSSG Rev 003 9/2008